

Code No.04-2901

## SPECIFICATION

## Dibasic Potassium Phosphate

Japanese Pharmaceutical Codex

REQUIREMENT	SPECIFICATION
* Description	White, crystals or mass
Identification	to pass test
pH	8.6~9.3
Purity	-
(1) Solubility	to pass test
(2) Chloride	max.0.009%
(3) Sulfate	max.0.017%
(4) Carbonate	to pass test
(5) Heavy metals	max.10ppm
(6) Sodium	to pass test
(7) Arsenic	max.2ppm
(8) Condensed phosphate	to pass test
Loss on drying	max.2.0%
Assay (after drying)	min.99.0%

Japanese Pharmaceutical Codex 2002(\*:Reference item)