

Code No.04-3235

## SPECIFICATION

N,N-Dimethylacetamide, Super Dehydrated  
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ nearly colorless, clear liquid
Solubility in ethanol	to pass test
Density (20°C)	0.938~0.943g/mL
Refractive index 20°C	1.435~1.441
Water	max.0.001%
Residue after evaporation	max.0.01%
Acidity (as CH <sub>3</sub> COOH)	max.0.005%
Assay (capillary GC)	min.99.5%