

## SPECIFICATION

## Diethyl Ether for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 300

REQUIREMENT	SPECIFICATION
* Appearance * Suitability for pesticide residue analysis (GC-FPD) Assay (C <sub>2</sub> H <sub>5</sub> OC <sub>2</sub> H <sub>5</sub> ) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis	Colorless clear liquid to pass test min.99.5%(area/area) 0.712~0.714g/mL max.0.1%(mass/mass) max.2ppm(mass/mass) to pass test

JIS K 8357:2023(\*:Additional test performed by Wako but not required by JIS)