

## SPECIFICATION

## Ethanol (99.5)

for Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.789~0.791g/mL
Refractive index 20°C	1.360~1.363
Absorbance (210nm)	max.0.70
Absorbance (220nm)	max.0.40
Absorbance (230nm)	max.0.20
Absorbance (240nm)	max.0.10
Absorbance (254nm)	max.0.02
Absorbance (260~400nm)	max.0.01
Water	max.0.1%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.9%