

SPECIFICATION

Ethyl Acetate for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
** Assay (CH ₃ COOC ₂ H ₅) (GC)	min.99.7%(area/area)
Density (20°C)	0.899~0.902g/mL
Water	max.0.05%(mass/mass)
Residue after evaporation	max.2ppm(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test

JIS K 8110:2023(*:Additional test performed by Wako but not required by JIS)(**:The specification is originally set by Wako not JIS)