

Code No.05-0701

SPECIFICATION

Ethanol(99.5) for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance * Suitability for pesticide residue analysis (GC-FPD) Assay (C ₂ H ₅ OH) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis	Colorless clear liquid to pass test min.99.5%(area/area) 0.789~0.791g/mL max.0.3%(mass/mass) max.2ppm(mass/mass) to pass test

JIS K 8093:2023(*:Additional test performed by Wako but not required by JIS)