

Code No.06-0419

SPECIFICATION

Formic Acid (abt.99%) for HPLC

| REQUIREMENT | SPECIFICATION |
|-----------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Solubility in water | to pass test |
| Absorbance (1→4, 254nm) | max.1.00 |
| Absorbance (1→4, 260nm) | max.0.20 |
| Absorbance (1→4, 270nm) | max.0.02 |
| Absorbance (1→4, 300~400nm) | max.0.01 |
| Water | max.1.0% |
| Fluorescence test | to pass test |
| Assay (HPLC) | min.99.5% |