

Code No.07-0589

SPECIFICATION

Geniposide

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
* Assay (qNMR) (JP method)	min.95.0%
** Solubility in methanol	to pass test
** Loss on drying over phosphorus(V) oxide in vacuum	max.2.0%

This product is suited for Geniposide for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia.
(*:The specification is originally set by Wako)(**):Additional test performed by Wako)