

Code No.08-0015

SPECIFICATION

Heptane

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|------------------------|
| * Appearance | Colorless clear liquid |
| Assay (C ₇ H ₁₆) (GC) | min.99.0%(area/area) |
| Density (20°C) | 0.682~0.686g/mL |
| Refractive index 20°C | 1.386~1.389 |
| Water | max.0.01%(mass/mass) |
| Residue after evaporation | max.0.001%(mass/mass) |
| Substances darkened by sulfuric acid | to pass test |

JIS K 9701:2024(*:Additional test performed by Wako but not required by JIS)