

SPECIFICATION

Hexane for Pesticide Residue And Polychlorinated Biphenyl Tests

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
Assay [CH ₃ (CH ₂) ₄ CH ₃](GC)	min.96.0%(area/area)
Density (20°C)	0.658~0.662g/mL
Water	max.0.01%(mass/mass)
** Residue after evaporation	max.1ppm(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test

JIS K 8825:2020(*:Additional test performed by Wako but not required by JIS)(**:The specification is originally set by Wako not JIS)