

## SPECIFICATION

Heptane, Super Dehydrated  
for Organic Synthesis

| REQUIREMENT                          | SPECIFICATION          |
|--------------------------------------|------------------------|
| Appearance                           | Colorless clear liquid |
| Density (20°C)                       | 0.682~0.686g/mL        |
| Refractive index 20°C                | 1.386~1.389            |
| Water                                | max.0.001%             |
| Residue after evaporation            | max.0.005%             |
| Substances darkened by sulfuric acid | to pass test           |
| Assay (capillary GC)                 | min.99.0%              |