

SPECIFICATION

Hexane, Deoxidized
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.001%
Residue after evaporation	max.0.001%
Acidity (as CH ₃ COOH)	max.0.002%
Dissolved oxygen	max.1 ppm
Fat	to pass test
Thiophene group	to pass test
Substances darkened by sulfuric acid	to pass test
Sulfide and reducing substances	to pass test
Assay (capillary GC)	min.96.0%