

Code No.08-1034

# SPECIFICATION

## Hyaluronic Acid Sodium Salt for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellow, crystals ~ powder or granules to pass test
Solubility in water	
Loss on drying at 60°C over phosphorus(V) oxide in vacuum	max.10.0%
Residue after ignition (as sulfate) (calculated on the dried basis)	15.0~20.0%
pH of a 1g/L solution at 25°C	6.0~7.0
Heavy metals (as Pb)	max.20ppm