

Code No.08-1039

SPECIFICATION

Heparin Sodium, from Porcine Intestinal Mucosa  
for Cell Culture

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale brown, powder or granules
Solubility in water	to pass test
Loss on drying at 60℃ over phosphorus(V) oxide in vacuum	max.10%
pH of a 10g/L solution at 25℃	5.0~8.0
Nitrogen compounds (as N) (after drying)	under investigation
Activity	min.130units/mg
Impurities	max.0.20%
Bacterial endotoxins	under investigation