

Code No.08-1049

## SPECIFICATION

## 4 Haloacetic Acid Mixture Standard Solution (t-Butyl Methyl Ether Solution)

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Impurities	to pass test
Chloroacetic acid (20°C)	950~1050mg/L
Dichloroacetic acid (20°C)	950~1050mg/L
Bromoacetic acid (20°C)	950~1050mg/L
Trichloroacetic acid (20°C)	950~1050mg/L