

## SPECIFICATION

## 2,3-Indolinedione

Guaranteed Reagent

| REQUIREMENT  | SPECIFICATION                    |
|--|----------------------------------|
| * Appearance   | Yellowish red, crystals ~ powder |
| Assay (C <sub>8</sub> H <sub>5</sub> NO <sub>2</sub> ) | min.97.5%(mass/mass)             |
| Solubility in N,N -dimethylformamide                   | to pass test                     |
| Melting point  | 200~204℃                         |
| Residue after ignition (as sulfate)                    | max.0.1%(mass/mass)              |
| Sensitivity  | to pass test                     |

JIS K 8089:2012(\*:Additional test performed by Wako but not required by JIS)