

Code No.09-0301

## SPECIFICATION

## 3-Aminomethyl-3,5,5-trimethylcyclohexylamine (mixture of cis- and trans-)

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Solubility in water	to pass test
Density (20°C)	0.920~0.928g/mL
Refractive index 20°C	1.485~1.491
Water	max.0.5%
Assay (capillary GC) (cis- + trans-)	min.98.0%