

Code No.09-0611

## SPECIFICATION

Interferon-  $\gamma$  , Human, recombinant, Animal-derived-free  
for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01ng/ $\mu$ g
Electrophoresis test	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source:E. coli expressed human interferon-  $\gamma$   
This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.