

Code No.11-0074

SPECIFICATION

Ketoprofen Standard for High Performance Liquid Chromatography

| REQUIREMENT | SPECIFICATION |
|----------------------------------|------------------------------|
| Appearance | White ~ nearly white, powder |
| Solubility in methanol | to pass test |
| Melting point | 94~97°C |
| Loss on drying at 60°C in vacuum | max.0.5% |
| Related substances | to pass test |
| Assay (HPLC) | min.98.0% |
| Assay (qNMR) | min.98.0% |