

## SPECIFICATION

## Lead(II) Oxide

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (PbO) Loss on ignition at 700°C Acetic acid-insoluble matter Chloride (Cl) Nitrate (NO <sub>3</sub> ) Sodium (Na) Potassium (K) Copper (Cu) Silver (Ag) Magnesium (Mg) Calcium (Ca) Bismuth (Bi) Iron (Fe)	Yellow ~ reddish yellow, powder min.99.5%(mass/mass) max.0.2%(mass/mass) max.0.02%(mass/mass) max.0.002%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.2ppm(mass/mass) max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass)

JIS K 8090 : 2019(\*:Additional test performed by Wako but not required by JIS)