

Code No.12-0104

## SPECIFICATION

### Linseed Oil

Wako 1st Grade

| REQUIREMENT           | SPECIFICATION                                                    |
|-----------------------|------------------------------------------------------------------|
| Appearance            | Slightly yellow ~ yellowish brown, clear ~ slightly muddy liquid |
| Density (20°C)        | 0.920~0.935g/ml                                                  |
| Refractive index 20°C | 1.480~1.484                                                      |
| Acid value            | max.1.0                                                          |
| Saponification value  | 188~196                                                          |
| Iodine value          | 170~204                                                          |