

## SPECIFICATION

L-Leucyl-p-nitroanilide Hydrochloride  
for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystalline powder ~ powder to pass test  +88~+92°  max.0.2%  to pass test  min.98.0%
Solubility in water	
Specific rotation [ $\alpha$ ]D20(c=1,H2O) (after drying)	
Loss on drying over sulfuric acid in vacuum	
Suitability for LAP measurement	
Assay	