

Code No.12-0409

SPECIFICATION

Lactitol Monohydrate

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|---|--------------------------------|
| Appearance | White ~ nearly white, crystals |
| Solubility in water | to pass test |
| Specific rotation [α]D20(c=10,H2O) | +13~+16° |
| Loss on drying at 90°C in vacuum | 4.5~6.5% |
| Residue after ignition (as sulfate) | max.0.1% |
| Assay (HPLC) | min.99.0% |