

Code No.12-0557

SPECIFICATION

L-Leucine Standard

for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in dil.hydrochloric acid	to pass test
Specific rotation [α]D20(c=4,HCl)	+14.7~+15.7°
Loss on drying at 105°C	max.0.2%
Suitability for amino acids determination	under investigation
Assay (volumetric analysis)	min.99.0%
Assay (subtracting method)	under investigation