

Code No.13-0193

SPECIFICATION

p-Methylaminophenol Sulfate
for Photography

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals
Solubility in water	to pass test
Solubility in hydrochloric acid	to pass test
Melting point (nitroso compound)	130~135°C
Residue after ignition (as sulfate)	max.0.1%
Diethyl ether extractive matter	max.0.2%
p-Aminophenol	max.2.5%
Heavy metals (as Pb)	max.0.002%
Iron (Fe)	max.0.005%
Assay	min.99.0%