

Code No.13-0747

SPECIFICATION

Methyl (R)-(-)-3-Hydroxybutyrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	1.056~1.061g/ml
Refractive index 20°C	1.420~1.424
Specific rotation [α]D20(neat)	-18.0~-23.0°
Water	max.0.2%
Optical purity	min.80.0%
Assay (capillary GC)	min.97.0%