

Code No.13-0884

## SPECIFICATION

Dichloromethane for Pesticide Residue And Polychlorinated Biphenyl Analysis  
Concentration 300

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH <sub>2</sub> Cl <sub>2</sub> ) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis * Suitability for pesticide residue analysis (GC- FPD) * Methanol (stabilizer)	Colorless clear liquid min.99.5%(mass/mass) 1.315~1.330g/mL max.0.01%(mass/mass) max.0.0002%(mass/mass) to pass test to pass test 0.2~0.5%(mass/mass)

JIS K 8117:2015(\*:Additional test performed by Wako but not required by JIS)