

# SPECIFICATION

## Magnolol

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

| REQUIREMENT                       | SPECIFICATION                        |
|-----------------------------------|--------------------------------------|
| Description (JP method)           | White crystals or crystalline powder |
| Identification (JP method)        | to pass test                         |
| Purity (JP method)                | -                                    |
| Related substances (JP method)    | to pass test                         |
| Unity of peak (JP method)         | to pass test                         |
| * Assay (qNMR) (JP method)        | min.95.0%                            |
| ** Solubility in methanol         | to pass test                         |
| ** Loss on drying over silica gel | max.1.0%                             |

This product is suited for Magnolol for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia.  
(\*:The specification is originally set by Wako)(\*\*):Additional test performed by Wako)