

## SPECIFICATION

## Magnoflorine Iodide

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay)

REQUIREMENT	SPECIFICATION
Description (JP method)	White to yellowish-pale gray, crystals or crystalline powder
Identification (1) (JP method)	to pass test
Identification (2) (JP method)	to pass test
Absorbance E1%1cm (223nm) (JP method)	1066~1132
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
* Assay (qNMR) (JP method)	min.95.0%
** Solubility in water-methanol	to pass test
** Water	measured value

\*:The specification is originally set by Wako

\*\*:Additional test performed by Wako