

## SPECIFICATION

Nitric Acid  
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (HNO <sub>3</sub> )	65~66%(mass/mass)
Residue after ignition (as sulfate)	max.3ppm(mass/mass)
Chloride (Cl)	max.0.3ppm(mass/mass)
Sulfate (SO <sub>4</sub> )	max.0.5ppm(mass/mass)
Copper (Cu)	max.0.1 ppm(mass/mass)
Lead (Pb)	max.0.1 ppm(mass/mass)
Arsenic (As)	max.0.01ppm(mass/mass)
Iron (Fe)	max.0.2ppm(mass/mass)

JIS K 8541:2021(\*:Additional test performed by Wako but not required by JIS)The density of this product is approximately 1.40g/mL