

Code No.14-0137

## SPECIFICATION

## Nitric Acid (1.38)

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Residue after ignition (as sulfate)	max.0.001%
Chloride (Cl)	max.1ppm
Sulfate (SO <sub>4</sub> )	max.5ppm
Copper (Cu)	max.0.5ppm
Lead (Pb)	max.0.5ppm
Arsenic (As)	max.0.05ppm
Iron (Fe)	max.1ppm
Assay	60~61%