

SPECIFICATION

Nickel(II) Nitrate Solution (10mg Ni/mL 15% Nitric Acid Solution)
for Atomic Absorption Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ green, clear liquid
Acidity (HNO ₃)	under investigation
Silver (Ag)	max.50ppb
Cadmium (Cd)	max.50ppb
Thallium (Tl)	max.100ppb
Lead (Pb)	max.100ppb
Arsenic (As)	max.100ppb
Selenium (Se)	max.100ppb
Assay (as Ni)	9.50~10.50mg/mL