

## SPECIFICATION

10 Nitrosamines Mixture Standard Solution (each 2  $\mu\text{g}/\text{mL}$  Methanol Solution)  
for Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Practical test	to pass test
N-Nitroso-N'-methylpiperazine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodimethylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosomorpholine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosomethylaminobutyric acid	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodiethylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosoethylisopropylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodiisopropylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosomethylphenylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodi-n-propylamine	1.97~2.03 $\mu\text{g}/\text{mL}$
N-Nitrosodi-n-butylamine	1.97~2.03 $\mu\text{g}/\text{mL}$