

SPECIFICATION

Oxolinic Acid Standard

for High Performance Liquid Chromatography

| REQUIREMENT | SPECIFICATION |
|--|---|
| Appearance | White ~ slightly yellow, crystals ~ powder to pass test |
| Solubility in sodium hydroxide solution | |
| Loss on drying over silica gel in vacuum | under investigation |
| Residue after ignition (as sulfate) | max.0.2% |
| Assay (HPLC) | min.99.0% |
| Assay (volumetric analysis) | min.99.0% |