

Code No.15-0277

SPECIFICATION

L(+)-Ornithine Monohydrochloride Standard for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystals ~ powder to pass test
Solubility in water	
Specific rotation [α]D20(c=5.5,H2O)	+10.0~+12.0°
Loss on drying at 105°C	under investigation
Suitability for amino acids determination	under investigation
Assay (volumetric analysis)	min.97.0%
Assay (subtracting method)	under investigation