

Code No.16-0005

## SPECIFICATION

## Palladium(II) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Brown ~ dark brown, powder
Assay (after drying)	min.99.0%
Solubility in dil. hydrochloric acid	to pass test
Loss on drying at 105 °C	max.1.0%
Nitrate	to pass test
Silver (Ag)	max.0.001%
Copper (Cu)	max.0.001%
Iron (Fe)	max.0.005%

JIS K 8154-1995

(\*:Additional test performed by Wako but not required by JIS)