

Code No.16-0005

SPECIFICATION

Palladium(II) Chloride

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--------------------------------------|----------------------------|
| * Appearance | Brown ~ dark brown, powder |
| Assay (after drying) | min.99.0% |
| Solubility in dil. hydrochloric acid | to pass test |
| Loss on drying at 105 °C | max.1.0% |
| Nitrate | to pass test |
| Silver (Ag) | max.0.001% |
| Copper (Cu) | max.0.001% |
| Iron (Fe) | max.0.005% |

JIS K 8154-1995

(*:Additional test performed by Wako but not required by JIS)