

Code No.16-0047

SPECIFICATION

Liquid Paraffin

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|--------------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 0.85~0.89g/ml |
| Acidity and alkalinity | to pass test |
| Chloride (Cl) | max.0.0015% |
| Sulfate (SO ₄) | max.0.002% |
| Sulfide | to pass test |
| Solid paraffin | to pass test |
| Substances darkened by sulfuric acid | to pass test |
| Substances reducing permanganate | to pass test |