

Code No.16-0102

## SPECIFICATION

## Phenol

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly red, crystals or mass
Assay (C <sub>6</sub> H <sub>6</sub> O)(GC)	min.99.0%(mass/mass)
Solubility in water	to pass test
Solubility in sodium hydroxide solution	to pass test
Freezing point	min.40.5°C
Residue after evaporation	max.0.02%(mass/mass)

JIS K 8798:2012

(\*:Additional test performed by Wako but not required by JIS)