

Code No.16-0107

## SPECIFICATION

## Phenolphthalein

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystalline powder
Assay (C <sub>20</sub> H <sub>14</sub> O <sub>4</sub> )	min.98.0%(mass/mass)
Solubility in ethanol	to pass test
Absorbance (0.01g/L, pH9.8)	min.0.60
Visual transition interval (pH)	7.8~10.0(colorless ~ red)

JIS K 8799:2020(\*:Additional test performed by Wako but not required by JIS)