

Code No.16-0112

SPECIFICATION

Phenol Red

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in ethanol Solubility in sodium hydroxide solution Absorbance (5mg/l, pH8.4) Loss on drying at 105°C Residue after ignition (as sulfate) Visual transition interval (pH)	Red ~ dark red, crystalline powder ~ powder to pass test to pass test min.0.53 max.3.0%(mass/mass) max.0.5%(mass/mass) 6.8~8.4(yellow ~ red)

JIS K 8800:2012(*:Additional test performed by Wako but not required by JIS)