

SPECIFICATION

Potassium Dichromate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (K ₂ Cr ₂ O ₇) Solubility in water Loss on drying at 105°C Chloride (Cl) Sulfate (SO ₄) Sodium (Na) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe)	Yellowish red ~ reddish yellow, crystals ~ crystalline powder min.99.5%(mass/mass) to pass test max.0.05%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.02%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.5ppm(mass/mass)

JIS K 8517:2006

(*:Additional test performed by Wako but not required by JIS)