

Code No.16-0382

## SPECIFICATION

## Potassium Hydrogen Phthalate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay [C <sub>6</sub> H <sub>4</sub> (COOK)(COOH)]	99.8~100.2%(mass/mass)
Solubility in water	to pass test
Loss on drying at 120°C	max.0.05%(mass/mass)
Chloride (Cl)	max.0.002%(mass/mass)
Sulfur compounds (as SO <sub>4</sub> )	max.0.006%(mass/mass)
Sodium (Na)	max.0.01%(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Lead (Pb)	max.5ppm(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8809:2020(\*:Additional test performed by Wako but not required by JIS)