

Code No.16-0429

## SPECIFICATION

## Dipotassium Hydrogenphosphate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ powder
Assay (K <sub>2</sub> HPO <sub>4</sub> ) (after drying at 110°C)	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	8.7~9.4
Loss on drying at 110°C	max.2.0%(mass/mass)
Chloride (Cl)	max.0.002%(mass/mass)
Nitrate	to pass test
Sulfate (SO <sub>4</sub> )	max.0.005%(mass/mass)
Sodium (Na)	max.0.5%(mass/mass)
Calcium (Ca)	max.0.03%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.001%(mass/mass)

JIS K 9017:2021(\*:Additional test performed by Wako but not required by JIS)