

Code No.16-0886

SPECIFICATION

Potassium Iodide for Oxidant Analysis
for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (after drying)	min.99.5%
Solubility in water	to pass test
Loss on drying at 110 °C	max.0.2%
Nitrogen compounds (as N)	max.0.001%
Heavy metals (as Pb)	max.5ppm
Suitability for oxidant determination	to pass test

JIS K 8827:2004
(*:Additional test performed by Wako but not required by JIS)