

Code No.16-1679

SPECIFICATION

Potassium Tetrafluoroborate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystalline powder ~
Solubility in water	powder
pH of a 50g/L solution at 25°C	to pass test
Chloride (Cl)	3.3~5.3
Sulfate (SO ₄)	max.0.01%
Heavy metals (as Pb)	max.0.01%
Iron (Fe)	max.0.001%
Assay	max.0.001%
	min.99.0%