

SPECIFICATION

Phenylfluorone

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in dil. hydrochloric acid-methanol Absorbance (7.5mg/L) Loss on drying at 105°C Residue after ignition (as sulfate) Suitability for tin determination	Yellowish red ~ brown, crystalline powder ~ powder to pass test min.0.70 max.5.0%(mass/mass) max.0.5%(mass/mass) to pass test

JIS K 9547 : 2018

(*:Additional test performed by Wako but not required by JIS)