

Code No.16-2121

## SPECIFICATION

### Polyoxyethylene(20) Sorbitan Monolaurate for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	Yellowish brown, clear ~ slightly muddy liquid
Viscosity (25°C)	300~500mPa·s
Solubility in water	to pass test
Water	max.5.0%
Acid value	max.3
Saponification value	40~50
Hydroxyl value	90~110
Deoxyribonuclease	to pass test
Ribonuclease	to pass test