

SPECIFICATION

Paeoniflorin

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Specific absorbance E1%1cm (230nm) (after drying)	245~285
* Water	measured value
* Assay (HPLC)	min.98.0%

*:Additional test performed by Wako