

## SPECIFICATION

## Polysorbate 20

## Japanese Pharmaceutical Excipients

REQUIREMENT	SPECIFICATION
* Description Identification Viscosity pH Sp. Gr. d 2020 Acid value Saponification value Purity (1) Heavy metals (2) Arsenic Water Residue on ignition	Slightly yellow ~ yellow, liquid to pass test 350~550mm <sup>2</sup> /s 4.0~7.0 1.090~1.130 not more than 4.0 43~55 - not more than 20ppm not more than 2ppm not more than 3.0% not more than 0.25%

Japanese Pharmaceutical Excipients(\*:Reference item)